

# PRESERVATION AND TRANSPORT MEDIA

## PreTect<sup>®</sup> TM

Product no: 400901



## DESCRIPTION

**PreTect TM** is a methanol based transport medium for preservation of mRNA in cytological samples, designed for use with the diagnostic test PreTect SEE and PreTect HPV-Proofer HV. The product enables preservation of mRNA for up to four weeks between 4°C to 30°C prior to isolation of nucleic acids.

**Packaging:** Box containing 96 vials (96x10 ml).

---

*Contact us for further information!*

---

### **PreTect AS**

Ustadhagan 8 | N-3490 Klokkearstua  
Norway | Tel: +47 32 79 88 00  
www.prectect.no | E-mail: prectect@prectect.no

 **PRETECT** | by |  **MEL-MONT MEDICAL, Inc.**

---

# Safety Data

According to 1907/2006/EC, Article 3

Date prepared: 28th March 2006

Revision: 10th February 2016

Version: 02-16/320101

## 1. Identification of the substance/mixture and of the company/undertaking

### 1.1 Product identifier

**Product Name:** PreTect TM  
**Article number:** 400101 / 400901/ 420101/420102

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses:** SU 22 Professional uses: Public domain (administration, education, entertainment, services, craftsmen)  
SU 20 Health services  
SU 24 Scientific research and development

**Uses advised against:** None

### 1.3 Details of the Supplier of the Safety Data sheet

**Manufacturer/Supplier:** PreTect AS  
Industriveien 8  
N-3490 Klokkestua  
Norway

**Telephone/Fax:** +47 32 79 88 00/+47 32 79 88 01  
**Email:** pretect@pretect.no

### 1.4 Emergency telephone number:

**Opening hours:** Only available during office hours

## 2. Hazards identification

### 2.1 Classification of the substance or mixture Classification according to 1272/2008 (CLP):

Flammable liquid (Flam. Liq. 3)  
Acute toxicity (Acute Tox. 4)  
Specific target organ toxicity - single exposure  
STOT SE 1

### 2.2 Label elements

**Hazard pictograms:**



**Signal word:**

Warning

Warning

Danger

### Hazard statements:

H226 Flammable liquid and vapour  
H302 Harmful if swallowed  
H312 Harmful in contact with skin

### Precautionary statements:

H332 | Harmful if inhaled.  
H370 | Causes damage to organs.  
P210 | Keep away from heat/sparks/open flames/hot surfaces. No smoking.  
P233 | Keep container tightly closed.  
P260 | Do not breathe dust/fume/gas/mist/vapours/spray.  
P280 | Wear protective gloves/protective clothing/eye protection/face protection.  
P301/  
P312 | IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.  
P307/  
P311 | IF exposed: Call a POISON CENTER or doctor/physician.