

# FOR ISOLATION OF NUCLEIC ACIDS

## PreTect® X

Product no: 111001 (96 isolations),  
111002 (960 isolations), 111004  
(bulk orders)



LYSIS

SiI

WB1

WB2

WB3

EluB



REF 111001

LOT X00XXYY

202Y-MM-DD



PreTect AS  
Ustadhagan 8  
Klokkestua  
NORWAY  
www.prelect.no

### PreTect AS

Ustadhagan 8 | N-3490 Klokkestua | Norway  
Tel: +47 32 79 88 00 | [www.prelect.no](http://www.prelect.no)  
E-mail: [prelect@prelect.no](mailto:prelect@prelect.no)



**PRETECT** | by | MEL-MONT MEDICAL, Inc.

## DESCRIPTION

**PreTect X** is a kit that uses the Boom method to isolate (purify and concentrate) total nucleic acids (RNA/DNA) from biological specimens prior to downstream amplification assay. (Validated for use with PreTect SEE and PreTect HPV-Proofer line products.)

### CE IVD

Includes Magnetic Silica, Wash Buffers 1, 2, 3, Elution Buffer. **Lysis buffer sold separately.**

N.B. PreTect X Wash Buffer 2 requires the addition of ethanol before use.